Passeo[®]-18 Lux[®]

<https://manuals.biotronik.com>

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Description

The Passeo®-18 Lux® paclitaxel-releasing PTA balloon catheter (hereinafter Passeo-18 Lux catheter, Passeo-18 Lux DCB, catheter or the device) is intended for dilatation of stenotic segments in infrainguinal arteries with simultaneous release of paclitaxel to the vessel wall in order to reduce the occurrence of a restenosis of the treated vessel segment. The balloon surface of the Passeo-18 Lux catheter is homogeneously coated with a delivery matrix incorporating 3 µg paclitaxel per mm² with a maximum amount of 15.2 mg on the largest balloon (7.0 x 200 mm). Paclitaxel is delivered to the vessel wall upon expansion of the balloon. The SafeGuard insertion aid protects the balloon in order to keep its factory-made profile and drug coating during insertion through the introducer sheath [see Figure 1]. The balloon is designed to inflate to a known diameter at a specific inflation pressure consistent with the compliance chart on the label. One radiopaque marker is located at each end of the balloon to facilitate angiographic visualization and positioning of the balloon catheter towards and across the lesion. The Passeo-18 Lux catheter includes a soft tapered tip to facilitate advancement of the device. The Passeo-18 Lux catheter has two luer ports at the proximal end. One port (inflation port) serves for connecting an inflation device to inflate/deflate the balloon. The other port enables flushing of the guide wire lumen. The Passeo-18 Lux catheter has a hydrophobic silicone coating on its outer surface. The Passeo-18 Lux catheter is compatible with guide wire and introducer sheath sizes according to the recommendations on the label.

How Supplied

Sterile. Non-pyrogenic. Device is sterilized with ethylene oxide.

Contents

- One Passeo-18 Lux catheter and one SafeGuard Insertion Aid with an Instruction Card in a sealed, peel-open pouch.
- One Instructions for Use manual.

Storage

Store protected from sunlight and keep dry. Store between 15°C [59°F] and 25°C [77°F]. Short-term storage temperature excursions between 10°C [50°F] and 40°C [104°F] are allowed up to 48 hours.

Indications

The Passeo-18 Lux catheter is indicated to dilate de novo or restenotic lesions in the infrainguinal arteries.

Contraindications

Contraindications for the device and peripheral dilatation catheters in general are:

- Lesions which cannot be reached or treated with the dilatation catheter.
- Allergy to contrast media, antiplatelet drugs or anticoagulants.
- Hemorrhagic diathesis or other conditions such as gastrointestinal ulceration or cerebral circulatory disorders restricting the use of platelet aggregation inhibitor therapy and anticoagulation therapy.
- Target stenosis located distal to a stenosis ≥ 50% that cannot be pretreated because the drug coating could get lost during crossing the proximal lesion.
- Allergy, intolerance or hypersensitivity to paclitaxel or structurally related compounds and/or to the delivery matrix n-butylryl tri-n-hexyl citrate (BTHC).
- Women who are breast-feeding, pregnant or are intending to become pregnant, or men intending to father children.

Warnings

- DO NOT use the device in the coronary, cervical and intracranial arteries.
- This device is designed and intended for single use only. DO NOT resterilize and/or reuse. Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient. Cleaning, disinfection, and sterilization may compromise essential material and design characteristics leading to device failure. The Manufacturer will not be responsible for any direct, incidental or consequential damages resulting from resterilization or reuse.
- DO NOT use the device if the outer package and/or inner package are/is damaged or opened, or if any of the provided information is obscured or damaged.
- DO NOT use the device after the Use by date indicated on the label.
- DO NOT expose the device to organic solvents e.g. alcohol.
- Use only an appropriate balloon inflation medium [e.g. 50:50 mixture by volume of contrast medium and saline]. NEVER use air or any gaseous medium to inflate the balloon.
- Only manipulate the device under angiographic guidance when it is inside the body.
- DO NOT advance or retract the device unless the balloon is fully deflated under vacuum. If strong resistance is felt during manipulation, stop the procedure and determine the cause before proceeding. Advancement by force may result in damage to the vessel and/or laceration or separation of the guide wire or dilatation catheter. This may necessitate recovery of device fragments.
- DO NOT exceed the rated burst pressure (RBP) indicated in the Compliance Chart. Use of a pressure-monitoring inflation device is mandatory to prevent over-pressurization.
- DO NOT exceed the original diameter of the vessel proximal and distal to the lesion when inflating the balloon to reduce the potential for vessel damage.
- DO NOT bend or squeeze the balloon portion of the device to prevent delamination of the drug coating.

Precautions

General Precautions

- Before using the device, the benefits and risks for each patient must be individually assessed.
- This device should only be used in medical facilities by physicians trained and experienced in performing vascular interventions (including cases of life-threatening complications).
- In order to prevent and reduce clotting, flush or rinse all devices entering the vascular system with sterile saline or a similar solution prior to use. The use of systemic heparinization during the procedure is recommended.
- After use dispose the device and packaging in accordance with hospital, administrative and/or local government policy.

Precautions Related to the Active Substance - Paclitaxel

- The amount of paclitaxel on the balloon surface corresponds approximately to a few tenths of the quantity usually used in antineoplastic treatment, which makes it rather improbable that interactions with other drugs will occur. However, caution should be exercised when concomitantly administering known CYP3A4 and/or CYP2C8 substrates (including terfenadine, cyclosporine, lovastatin, midazolam, ondansetron) or drugs with high PPB (especially sulfonyleureas, anticoagulants of the coumarin type, salicylic acid, sulfonamides, digitoxin). For possible interactions with other drugs within the scope of paclitaxel administered for oncological indications, the relevant instructions for use should be consulted. A study of possible interactions with other drugs within the scope of paclitaxel in association with treatment accompanying drugs has not been established.
- The implantation of a drug-eluting stent or another drug releasing balloon catheter at the same site should be avoided since an overdosage or interaction with the active agents cannot be excluded. The use of devices as indicated above in a new subsequent intervention respectively has not been evaluated.
- In case of the treatment of long lesions (longer than max. balloon length available) the lesion can be treated with multiple Passeo-18 Lux catheters. Each additional Passeo-18 Lux catheter should be used with 10 mm overlap.
- It is recommended to use gloves, mouth, nose and eye protection in the unlikely event that the active compound of the balloon coating is released prior to insertion of the Passeo-18 Lux catheter into the introducer sheath. Potentially the released particles may enter the respiratory tract.
- The treating physician should balance the medical advantage of treating a patient with a paclitaxel-releasing balloon against the risk of potential adverse events as described in section "Potential Adverse Events/Complications".
- Paclitaxel-related systemic adverse events as described in the "Potential Adverse Events/Complications" section are not expected after treatment with a single Passeo-18 Lux catheter. If more than one catheter is used in the same procedure, consider that the safety of a total drug dosage exceeding 31.2 mg in a patient has not been clinically evaluated. See table "Sizes and Associated Maximum Drug Loads" for reference.
- Possible exposure of the user to paclitaxel can be minimized when the Passeo-18 Lux catheter is used according to the IFU. Paclitaxel is bound to the expient BTHC on the balloon surface and it is unlikely that it propagates through the air (e.g. like a powder or dust) and beyond a radius of 1 m. Nevertheless, all people handling the Passeo-18 Lux catheter should wear gloves, mouth, nose and eye protection and avoid any unprotected direct or indirect contact (e.g. through contaminated material or liquids).
- Hospitals and personnel handling the Passeo-18 Lux catheter should take all of the appropriate precautions for handling and disposal of the product (including packaging) and be trained in the handling of cytostatic drugs.
- Paclitaxel is in potential a genotoxic (in particular aneugenic/clastogenic) agent based upon its pharmacodynamic mechanism of action, which is interference with the microtubule disassembly. The relevance of this specific mechanism of genotoxicity for human carcinogenicity risk is currently not known.

Precautions Related to Special Patient Populations and Indications

- The risk associated with antiplatelet treatment should be taken into account. Special consideration is required in patients with recent active gastritis or peptic ulcer disease (PUD).
- It is not possible to state the effects of the Passeo-18 Lux catheter on the unborn child. No clinical data is available on the use of the Passeo-18 Lux catheter in pregnant women or men intending to father children, i.e. reproduction relevant contraindications and risks are unknown.
- Pediatric use: The safety and effectiveness of the Passeo-18 Lux catheter in pediatric patients have not been established.
- Data from literature (human and animal) has shown transfer of paclitaxel to breast milk and a potential toxic effect in a nursing infant cannot be excluded.

Procedure Precautions

- Prior to the procedure, the device should be visually examined to verify its integrity and to ensure that its size is suitable for the specific procedure for which it is to be used.
- DO NOT use the device if the Safeguard insertion aid does not fully cover the balloon after removal from the protection ring.
- It is recommended to use gloves, mouth, nose and eye protection in the unlikely event that the active compound of the balloon coating is released upon removal of the device from the protection ring. Potentially, the released particles may enter the respiratory tract.
- Pre-treatment of significant stenosis proximal to target lesion is required to prevent delamination of the balloon coating during crossing the lesion.
- Exercise care when handling the device to reduce the possibility of accidental breakage, bending or kinking of the catheter shaft.
- Use only guide wires with a maximum 0.018" (0.46 mm) diameter.

- Use only the minimum introducer sheath sizes stated on the label. However, if Passeo-18 Lux catheter is used in conjunction with long and/or braided introducer sheaths, a larger French size may be necessary to reduce friction.
- Contact with liquids prior to insertion, touching of the balloon or wiping of balloon surface should be strictly avoided, since this could wash away or damage the balloon coating.
- The inflation diameter of the balloon must not exceed the diameter of the artery proximal or distal to the stenosis.
- Keep the SafeGuard insertion aid in place while preparing and purging air from the device.
- The raised portion of the SafeGuard insertion aid must remain partially visible outside the introducer during insertion.
- Once the proximal end of the balloon is within the introducer, slide back the SafeGuard insertion aid out of the introducer sheath and towards the hub of the device.
- DO NOT use the device if vacuum cannot be held as this indicates a leak in the system.
- Balloons exceeding 3.5 mm diameter and 60 mm length should be deflated for at least 90 seconds to ensure that all inflation medium was completely removed.
- It is recommended to maintain negative pressure whenever the device is withdrawn.

Potential Adverse Events/Complications

Potential adverse events associated with the included drug (paclitaxel) and delivery matrix (BTHC):

- Allergic/immunologic reaction to the drug (paclitaxel or structurally related compounds) or to the balloon delivery matrix n-butylryl tri-n-hexyl citrate (BTHC).
- Alopecia.
- Anemia.
- Gastrointestinal symptoms.
- Hematopoietic disorders [incl. leukocytopenia, neutropenia, thrombocytopenia].
- Hepatic enzyme changes.
- Histological changes in the vascular wall, incl. inflammation, cell damage or necrosis.
- Disorders of the heart conduction system.
- Myalgia / arthralgia.
- Peripheral neuropathy.
- Pseudomembranous colitis.
- Mutagenic effects.

Note: The very small quantity of paclitaxel in the blood plasma means that the classic undesirable effects caused by paclitaxel appear less relevant than for a systemic treatment. However, side effects yet unknown cannot be ruled out.

Potential adverse events associated with PTA procedure include, but are not limited to:

- Balloon catheter events: failure to reach or cross the lesion, inflation difficulties, deflation difficulties, withdrawal difficulties.
- Vascular events: pseudoaneurysm, arteriovenous fistula formation, vessel dissection, rupture or perforation, vessel recoil or restenosis, thrombosis or occlusion, vasospasm, peripheral ischemia, embolization, local inflammatory reactions resulting from possible injuries of the vascular wall.
- Neurologic events: peripheral nerve injury / neuropathy.
- Bleeding events: access site hematoma, access site bleeding, hemorrhage requiring transfusion or other treatment.
- Accompanying medication related events: side effects per corresponding package leaflet.
- Allergic reactions to contrast media, antiplatelet drugs or anticoagulants.
- Infection.
- Tissue necrosis and limb loss.
- Death.

Furthermore, all procedure-related adverse events as described in the national and international guidelines of the respective medical associations apply.

Medication Schedule

Appropriate anticoagulant, antiplatelet and vasodilation therapy should be administered to the patient.

The following medication regimen is solely intended as a possible guideline and should not be seen as a strict guidance.

Pre-procedure medication regimen:

- ASA 300 mg, 12 hours before procedure.
- Clopidogrel 300 mg, 12 hours before procedure

Procedural medication regimen:

- Intraarterial bolus of heparin (3000 - 5000 U).

Post-procedure medication regimen:

- ASA 100 mg daily, indefinitely.
- Clopidogrel 75 mg daily, for 4 weeks (check leukocyte count and platelet count after 2 and 4 weeks).

Directions for Use Device Selection

Select an appropriate balloon for the target vessel. The diameter of the balloon should closely match the reference diameter of the target vessel. Select a balloon length which closely matches the lesion length. If the lesion cannot be crossed with the desired Passeo-18 Lux catheter, a smaller diameter non drug-coated PTA catheter must be used or appropriate vessel preparation performed before finally dilating with a Passeo-18 Lux catheter.

In case of the treatment of long lesions (longer than max. balloon length available) the lesion can be treated with multiple Passeo-18 Lux catheters. Each additional device should be used [sequentially] with 10 mm balloon overlap [see Figure 2].

Caution: For prevention of local overdosing it is not indicated to use a second Passeo-18 Lux catheter or any other drug-coated balloon at the same treatment site, respectively. Also the implantation of a drug eluting stent at the same site should be avoided since an over-dosage or interaction between the active agents cannot be excluded.

Device Preparation

1. Take the protection ring with the device out of the package and place it onto a sterile field.
2. Gently pull out the catheter from the protection ring.

Caution: DO NOT remove the SafeGuard insertion aid yet, in order to avoid any contact with the balloon coating.

Flushing of the Guide Wire Lumen

3. Connect a 10 ml or 20 ml syringe containing sterile saline to the luer port of the guide wire lumen at the proximal end of the device.
4. Flush the guide wire lumen.
5. Remove the syringe.

Removal of Air from the Inflation Lumen of the Device with the Insertion Aid still in place

6. Fill a 20 ml capacity inflation device with 6 ml of inflation medium [e. g. 50:50 mixture by volume of contrast medium and saline]. Use 10 ml of inflation medium for balloons exceeding 5 mm in diameter and 120 mm in length.
7. Remove air from the inflation device according to the manufacturer's recommendations and instructions.
8. Attach the inflation device to the inflation port of the catheter.
9. Pull back on the plunger and aspirate for 30 seconds until no bubbles appear inside the barrel of the inflation device. Repeat this process several times if necessary.
10. Return to neutral pressure.

Preparation of SafeGuard Insertion Aid

1. Advance the balloon inside the SafeGuard insertion aid to ensure the distal tip of the balloon is visible beyond the distal end of the SafeGuard insertion aid [see Figure 3]. DO NOT REMOVE the SafeGuard insertion aid.

Insertion Technique

2. If not already done, position the guide wire, under angiographic guidance, in accordance with standard PTA techniques.
3. Thread the distal tip of the catheter over the proximal end of the guide wire and advance until the guide wire exits the luer lock at the proximal end of the catheter.
4. Holding the device by the SafeGuard insertion aid [see Figure 4], carefully insert the Passeo-18 Lux catheter together with the SafeGuard insertion aid into the introducer sheath, until the raised portion of the SafeGuard insertion aid contacts the valve of the introducer sheath [see Figure 5].
5. With one hand, continue holding the SafeGuard insertion aid to maintain its position. With the other hand, advance the device until the proximal end of the balloon is inside the introducer sheath.
6. While holding the device in place, slide back the SafeGuard insertion aid out of the introducer sheath and towards the hub of the catheter [see Figure 6].
7. If necessary, the SafeGuard insertion aid can be peeled off by pushing it onto the kink protector and then gently pulling it away from the kink protector at the proximal end of the shaft [see Figure 7].
8. Advance the device over the guide wire towards the lesion.
9. Position the balloon across the lesion using the balloon radiopaque markers as reference points.

Balloon Inflation

20. Inflate the balloon with the inflation device according to the Compliance Chart to dilate the lesion using standard PTA techniques for at least 30 seconds.
21. If a significant stenosis persists, inflate the balloon again increasing pressure gradually until the lesion fails to improve.

Note: Administration of adequate drug contents occurs only during the first inflation.

Balloon Deflation

22. To deflate the balloon completely, pull back the plunger of the inflation device and lock it in this position. Apply vacuum to the balloon under angiographic control for at least 30-90 seconds depending on the balloon size.

Device Withdrawal

23. While maintaining vacuum with the inflation device and stable guide wire position across the lesion, withdraw the catheter carefully from the lesion and out through the introducer sheath. If the balloon has been dilated several times, there may be some resistance as it is pulled back into the introducer sheath. In case of difficulties, pull out the catheter and the introducer sheath together.
24. Completely remove the catheter from the guide wire.
25. After use, the device and packaging may be contaminated with infectious substances (e.g. blood). Therefore, dispose the device and packaging in accordance with hospital, administrative and/or local government policy.

Warranty/Liability

This device and each of its components (hereinafter device) were designed, manufactured, tested and packaged with all reasonable care. However, since the Manufacturer does not have control over the conditions under which the device is used, the contents of this Instructions for Use (IFU) are

to be considered as an integral part of this disclaimer for cases when a disturbance of the device's intended function may occur for various reasons.

The Manufacturer does not guarantee that the following events will not occur:

- Device malfunctions or failures
 - Patient's immune response to the device
 - Medical complications during the use of the device or as a consequence of the device being in contact with the patient's body
- The Manufacturer does not assume any liability for:
- The use of the device that is not in accordance with the stated intended use/ indication, contraindications, warnings, precautions and the directions for use of this IFU
 - Modification to the original device
 - Events which could not have been foreseen at the time of device delivery using the available levels of science and technology
 - Events originating from other Manufacturer's devices or devices from other manufacturers
 - Force majeure events including, without being limited to, natural disasters.
- The above provisions shall be without prejudice to any disclaimer and/or limitation of liability agreed separately with the customer to the extent permitted by applicable laws.

Intellectual Property

The following list includes trademarks or registered trademarks of the BIOTRONIK Group of Companies in the EU, the United States, and possibly in other countries: BIOTRONIK, Passeo, Lux, SafeGuard. All other trademarks are the property of their respective owners.

Safety and Clinical Performance

The Summary of Safety and Clinical Performance (SSCP) can be found at <https://ec.europa.eu/tools/eudamed> using the Basic UDI-DI: 76401304BUDI00008GK

Clinical Benefits

Dilation of stenotic segments in infrainguinal arteries with simultaneous release of paclitaxel to the vessel wall is expected to result in reduction of repeat revascularization and/or preservation of a functional limb.

Reporting of Serious Incidents

Report any serious incident that has occurred with the device to the Manufacturer and to the competent authority of the country in which your hospital is established. In case of a device malfunction, stop using the device immediately and return it to the Manufacturer.

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All rights reserved. Specifications subject to
© 2025 BIOTRONIK AG
BIOTRONIK AG